

This certificate is an example. A complete certificate containing lot-specific information, such as concentration and volume, will be shipped with your purchase.



PAN-Chromeo™ 546 Proteinic-blocked (PAN Nanoparticles)

Catalog No: 15231, 16231

Lot No: Not Available

Chemical Properties:

Contents: 5 ml (Cat # 15231) or 10 ml (Cat # 16231) aqueous suspension of PAN-Chromeo 546 Proteinic blocked, 0.5 % (w/w) in 10 mM MOPS + 0.5 % (w/w) BSA, pH 7, 15mM sodium azide.

Fluorescent Properties:

PAN-Particles are based on a novel polymer, which form lattices with a reduced oxygen permeability and results in a very good protection of the incorporated dye against other environmental effects. The combination of this new polymer with the Chromeo-dyes results in novel nanoparticles with high fluorescent intensities. The particles are smaller than 50 nm, the main fraction is about 30 nm in diameter. The surface is carboxy-modified, has additional sulfo-groups and is blocked with BSA. PAN-Nanoparticles are "loaded" with the lipophilic-form of Chromeo 546, showing similar fluorescent properties as the yellow-orange Cy3 dye. Filters designed for Cy3 can be used with Chromeo 546 carboxylated particles.

Excitation Max: 520, 540 nm

Emission Max: 577 nm

Diameter (average): 30 nm

Delivery-Form: 0.5 % (w/w) in buffer

Quality Control:

The particles have been quality tested by spectro-photometrical evaluation and by determining the solid residue after lyophilization.

Storage:

To ensure stability, the aqueous suspension should be stored at 4°C in the dark.

Caution: Do not dry or freeze; do not use more than 5% organic solvent. This product is guaranteed for 12 months from the date of arrival.

